


# SECUTAPE®

## FIXATIONS Since 1986



# CATALOG



**TechniMed AG**

SWITZERLAND since 1987

fon 0041 / 71 844 49 20

Kirchstrasse 54

CH-9400 Rorschach

fax 0041 / 71 845 52 45

**Germany**

fon 0800 007 9892

fax 0800 941 3012

**Austria**

fon 01/ 440 7979

fax 01/ 440 7990

mail: [info@technimed.ch](mailto:info@technimed.ch)

**CE**

Made in Germany

[www.secutape.ch](http://www.secutape.ch)

Copyright by TechniMed AG


# SECUTAPE®

**For 40 years now, SECUTAPE is offering:**

field-proven, high-quality, skin-friendly, and very secure fixations for **all tubes, sonds and lines on the human body.**

Our idea of sticking a plaster to the skin and attaching a tab made of Velcro and adhesive to it was registered as our first patent in 1986.

**Back in 1986, we launched the world's first mass-produced fixations.**



**Our standard fixations 2 and 3.**

They have proven themselves millions of times over and have become the standard for many applications in clinics.

Gradually, we have patented further innovative solutions and brought them to market. Today, we offer **special fixatons for all connections and lines on t the body.**

We manufacture our products in ergonomic shapes from modern high-tech materials **in Germany and Switzerland.**

Thanks to the use of innovative, high quality materials developed specifically for us, SECUTAPE adheres better and longer, which also makes it **cost-effective over time.**

Of course, we are copied, but we are somehow proud with the fact that only good products are copied.

Our Fixations are moustly single packaged after the EU MDR regulations. They all have an EAN code as well as an QR code that links to an application video.

**For the EU we bill and deliver out of Germany**

# SECUTAPE®

## Size 2+ / 3+ Standard fixation


Since over 35 years our reliable fixations!

for the fixation of :

- PEG
- CAPD
- CVC
- suprapubic catheters
- infusions
- urin-catheter
- redondrains
- ports
- wounddrains
- and all other tubes

For heavy tubes pls. use size 3+


The fixation can be opened and closed many times.




 **ord.no. 8500201** 1 PU 200 pcs.  
NW basic  
singlepack 8 bags a.25pcs.

**2+** **ord.no. 8500211** 1 PU 1000 pcs.  
20 bags. à 50 pcs.  
9X3.5 cm




 **ord.no. 8500301** 1 PU 200 pcs.  
NW basic  
singlepack 8 bags a.25pcs.

**3+** **ord.no. 8500311** 1 PU 1000 pcs.  
20 bags. à 50 pcs.  
15x4.5cm



Stays appr. 3 - 5 days



# SECUTAPE®

## Velcro-binder

### A new universal system for all applications and diameters!

Velcro binder holds reliably, **without any adhesive** tubes and wires with the new patented NANOPLAST® antislip foil.

The fixation can be opened and closed many times without leaving any marks on the tube!


Velcro binder comes as a set with different holders, also in Hydrocolloid for every application.

Due to the high quality breathing material, the fixation stays even when it gets wet.


NANOPLAST® has the so-called Gecko effect and is made out of silicon with aprox. 29000 nano legs per cm .

The fixation is very good for urine catheter and especially for the fixation of implanted tubes/lines.


### Velcro-binder



ord.no. 8501000 PU 200 pcs.  
singlepack. 8 bags a 25 pcs.




Velcro binder - Set SMALL  
**Hydrocolloid**  
ord.no. 8501100 PU 100 pcs.  
singlepack. 4 bags a 25 pcs.



Velcro binder - Set BIG  
ord.no. 8501130 PU 100 pcs.  
singlepack 4 bags a 25 pcs.

**Stays up to 3 days,  
Hydrocolloid up to 5 days**



CE

Made in Germany


[www.secutape.ch](http://www.secutape.ch)

## Fixation for "tubes and wires" especially for CVC

UNIFIX is perfectly suited for the fixation of thinner catheters, Probes and lines, e.g. multi-lumen CVC...

The basis is a very good holding, ergonomic shaped holder made of skin-friendly, breathing tape that also holds when it gets wet.

The catheter is placed in an adhesive bed and held securely with a Velcro strap.



UNIFIX can also be used to attach filter, transducer or similar parts securely to the body.



### Unifix

ord.no. 8501500 PU 200 pcs.

8 bags. à 25 pcs.  
singlepack

holder 6 x 6 cm



### Unifix Hydrocolloid

ord. no. 8501509 PU 200 pcs.

8 bags. à 25 pcs.  
singlepack

holder 6.5 x 5 cm




As an additional fixation for the tubes a non-slip holder with a clamp.

### KROKOFIX small

ord. no. 8501410



The use of two fixations reduces strain and movement on the CVC, which prevents infections.





## Fixing set for big lumina

### Our best holding Fixation for large tubes!


Good for mobile patients!

#### Specially designed for:

chest tubes, ECMO, rinsing catheters and  
for the fixations of urinary catheters.

Through a special coating it even holds  
silicon catheters!

Due to the high quality nonwoven material,  
**the fixation stays even when it gets wet.**



**ord.no. 8501200 1 PU 100 pcs.**  
**hydrocolloid** 4 bags. à 25 pcs.  
singlepack.  
12x6.5cm



**ord.no. 8501201 1 PU 100 pcs.**  
**nonwoven** 4 bags. à 25 pcs.  
singlepack.  
12x6.5cm



**Stays appr. up to 4 days  
hydrocolloid up to 6 days.**



## Nasal fixations

Fixation for nasale sondes

### Patented two piece fixation for all nasal sondes.

A specially shaped basic holder made of a skin friendly and breathable nonwoven or a hydrocolloid.

The basic holder has a Velcro part, a second strap with Velcro is glued to the tube.

The holder remains on the nose, if necessary only the straps are changed.

**There is no more skin damage** caused by the constant tearing of the plaster strips. The probe is hanging loose in the nose without causing any problem.




**ord.no.. 8500400** **1 PU 200 pcs.**  
 hydrocolloid 8 bags a.25pcs.  
 singlepack.  
 set  
 holder 4.4 x 3.5cm  
 strap 7.5 x 0.8 cm



**ord.no.. 8500401** **1 PU 200 pcs.**  
 nonwoven 8 bags a.25pcs.  
 singlepack.  
 set  
 holder 4.4 x 3.5cm  
 strap 7.5 x 0.8 cm



**ord.no.. 8500410** **1 PU 400 pcs.**  
 board à 4 pcs.  
 singlepack.



**ord.no.. 8500440** **1 PU 500 pcs.**  
 board à 5 pcs.  
 singlepack.

**Stays appr. 3 days  
hydrocolloid 5 days**

# SECUTAPE®

## Fixing set for endotracheal tubes oral/nasal

**Since over 30 years our reliable fixation for endotracheal tubes!**

Easy to fit by just one person and inexpensive due to longer application time..

Two holders of a very good adhering and skin friendly nonwoven or hydrocolloid with Velcro

Non stretchable, hypoallergenic nylon straps which do not absorb liquid.  
The straps are centrally coated with a special adhesive.

**The system holds very well when used correctly and allows safe aspiration without a second person.**

The standard strap is **soft on one side.**

The Softband is **soft on both sides**, which prevents skin damages.



**Stays appr. up to 3 days,  
hydrocolloid up to 5 days.**



Made in Germany  
[www.secutable.ch](http://www.secutable.ch)

# SECUTAPE®

## Tube fixing set EXTRA


**Our known tube fixing set with a newly developed strap.**

With our proven nonwoven holders and a very soft hypoallergenic nylon strap that does not absorb humidity and **does not stretch**.

The band has a slot in the center for a secure position of the tube.


The band is very well held on the tube with a self-adhesive Velcro strip.

The system holds very well when correctly installed and allows safe aspiration **without a second person**.



**Fixing set EXTRA**  
for endotracheal tubes

**ord.no. 8500610** 1 PU 200 pcs.  
singlepack 8 bags à 25 pcs.  
Both parts are available separately



CE

Made in Germany

[www.secutape.ch](http://www.secutape.ch)



## Short-term fixing tape

self-adhesive  
for endotrachealtubes.

### Fixing tape for the short-term fixation of endotracheal tubes.

Newly developed high-quality self-adhesive fixing tape,  
specially designed for the anesthesia.

The hypoallergenic, skin-friendly and excellent  
holding tape does not stretch and rip.

It is packed in a hygienic dispenser and easy to handle  
due to its shape and the special liner.

The perforated dispenser is opened at the front, which  
eases the withdrawal and storage.



### Short-term fixing tape

for endotracheal tubes

**ord.no. 8500820** 1 PU 400 pcs.  
box a. 100 pcs.  
37x3cm



This strap does not replace  
a good long-term fixation!



CE

Made in Switzerland  
[www.secutape.ch](http://www.secutape.ch)

# SECUTAPE®

## Tube-band


### Pat. fixation for endotracheal tubes

A hypoallergenic soft nylon band, **it does not stretch** and does not absorb any liquid.

Since the strap does not stretch no need to tight.

Due to the silicon coating on one side the strap has no hard edges, does not slip on the skin and the tube stays in place!

It is easy in handling, prevents most skin problems and must be changed less frequently than conventional textile or foam straps.



**Stays approx. up to 2 days.**



CE

Made in Germany  
[www.secutape.ch](http://www.secutape.ch)

# SECUTAPE®

## Universal-band

### For endotracheal tubes & tracheal cannulas

Newly developed double-sided nylon band is hypoallergenic and pleasantly soft.

**It does not stretch** and does not absorb liquid. Therefore the patient does not have a humid band around the neck.

The band is also easy to clean.

It has to be changed less frequently than the usual foam or fabric bands.

Since the band should not be applied tight, there will be less skin damage and more comfort for the patient.

**One product for both applications!**



### Universal-band

ord.no.8500890 1 PU 200 pcs.  
singlepack. 8 bags à 25 pcs.

length 75 cm



Time of use approx. up to 3 days.



CE

Made in Germany

[www.secutape.ch](http://www.secutape.ch)



**Tracheofix**  
Comfort


## The advantage of new materials!

### Fixation band for tracheal cannulas

Newly developed thick and fluffy nylon band.  
Hypoallergenic and pleasantly soft.

**The high quality band does not stretch and does not absorb liquid.**

It is easy to apply and due to its material it has to be changed less frequently than the usual foam or fabric bands.



### Superabsorber Pad

In the case of skin damage and heavily sweating patients a superabsorber pad is available.

The med. grade absorber pad sticks by itself and keeps the skin dry.

It is therefore an ideal solution for long-term fixation.



### Tracheofix Comfort

length 50cm x 3cm

ord.no. 8501750 PU 100 pcs.

singlepack. 4 bags. a 25 pcs.



### Superabsorber Pad

length 20 cm

ord.no. 8501770 PU 100 pcs.

singlepack. 4 bags. a 25 pcs.



**Stays up to 3 days  
including absorber pad  
up to 5 days.**



Made in Germany  
[www.secutable.ch](http://www.secutable.ch)



# Krokofix

## Fixation for "tubes and wires" at clothes and bed linen.

- two functions!
- simple handling !
- two sizes with stainless steel clamp.
- secure and non slipping adjustment!
- variable use!
- disinfection with usual means!

A manual anti-slip fixation for tubes of all sizes.  
Good for mobile patients and strong pulling.



**Just clamp it to clothes or bed linen.**



**Krokofix big**  
ord.no. 8501420 1 PU 200 pcs.  
23 cm



**Krokofix small**  
ord.no. 8501410 1 PU 200 pcs.  
16.5cm



**SECUTAPE®**

## Head holder-set ECMO

### Adjustable head holder set for fixation of an veno-arterial ECMO

The elastic head holder is a thick and soft nylon band and fits all head sizes.

The holder is, in combination with the Velcro straps and the Hydrocolloid holders, a very reliable Fixation, which gives the patient a certain room for movement.


The additional Fixation on the neck helps to reduce the strain on the hoses.

With the enclosed KROKOFIX, which should be fixed to the bed linen, to reduce the strain on the drains.



### Head holder-set: ECMO

ord.no. 8501610, PU 10 set



CE


Made in Germany

[www.secutape.ch](http://www.secutape.ch)

**SECUTAPE®**

**FIXIERUNG  
FIXATION  
FISSAGGIO**

*SIMPLY SAFE*



[www.secutable.ch](http://www.secutable.ch)